# SIDICOM® WLS V 3.0

Konfigurationsanleitung / Bedienungsanweisung  
Configuration Instructions / Operating Instructions

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# SIDICOM® WLS V 3.0

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The screenshot shows the SIDEXIS DICOM WLS application window. The top menu bar includes 'Query', 'Order', 'View', and 'Help'. Below the menu is a toolbar with icons for patient selection, search, and navigation. The main window is divided into two sections. The upper section displays a table of patient worklist entries. The lower section shows the detailed DICOM metadata for the selected patient, Ballack^Michael.

Order	Order Time	Patient (Name, date of birth)	Patient's ID	Modality	Image Type	Accession no.
1	20.07.07 12:00:00	Friedrich^Arne 29.05.79	09297264	OT		99750715
2	20.07.07 12:00:00	Ballack^Michael 26.09.76	23793809	IO		02003845
3	20.07.07 12:00:00	Lehmann^Jens 10.11.69	31224737	CT		00466792
4	20.07.07 12:00:00	Hildebrand^Timo 05.04.79	01288080	IO		51822021
5	20.07.07 12:00:00	MeierHuber^Fred 10.11.90	47550557	IO	IO	93306300
6	20.07.07 12:00:00	Hans^Peter	49454944	OT		25040326
7	20.07.07 12:00:00	Meier^Hans 04.03.87	50209415	CT	3D	49849553
8	20.07.07 12:00:00	Lukas^Podolski 04.06.85	20678302	DX		69596292
9	20.07.07 12:00:00	Kahn^Oliver 15.06.69	15006754	CT		13546521
10	20.07.07 12:00:00	Klose^Miroslav 09.06.78	48113680	DX		04113786

ID:(0008,0000)	IdentifyingGroupLength =	40
ID:(0008,0050)	AccessionNumber =	02003845
ID:(0008,0080)	InstitutionName =	
ID:(0008,0090)	ReferringPhysiciansName =	
ID:(0008,1110)	ReferencedStudySequence :	
ID:(0010,0000)	PatientGroupLength =	152
ID:(0010,0010)	PatientsName =	Ballack^Michael
ID:(0010,0020)	PatientID =	23793809
ID:(0010,0030)	PatientsBirthDate =	19760926
ID:(0010,0040)	PatientsSex =	
ID:(0010,1000)	OtherPatientIDs =	
ID:(0010,1001)	OtherPatientNames =	
ID:(0010,1020)	PatientsSize =	1.89
ID:(0010,1030)	PatientsWeight =	84.0
ID:(0010,2000)	MedicalAlerts =	
ID:(0010,2110)	ContrastAllergies =	
ID:(0010,2154)	PatientsTelephoneNumbers =	
ID:(0010,2160)	EthnicGroup =	
ID:(0010,21b0)	AdditionalPatientHistory =	
ID:(0010,4000)	PatientComments =	
ID:(0020,0000)	ImageGroupLength =	64
ID:(0020,000d)	StudyInstanceUID =	1.2.40.0.13.0.192.168.15.58.18481629.1147082817253.32772

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Bestell-Nr. **62 01 250 D 3348**

# SIDICOM® WLS V 3.0

Configuration Instructions/ Operating Instructions

English

### Contents

I	Overview .....	4
I.1	Structure of these instructions .....	4
I.2	Symbols used in this manual .....	5
I.3	Abbreviations and Definitions .....	6
I.4	Focus .....	7
I.5	System Requirements .....	7
I.6	Scope of delivery of SIDICOM WLS .....	7
I.7	New features of the current version .....	8
II	Operating Instructions .....	10
II.1	Workflow – Order based exposure .....	10
II.1.1	Requesting order data and patient demographics .....	10
II.1.2	Selecting and activating an order .....	13
II.1.3	Executing an order .....	16
II.2	Workflow – Creating a new DICOM study .....	20
II.3	Link tests .....	21
II.4	Send tests .....	21
III	Installation and configuration instructions .....	23
III.1	What is required prior to installing SIDICOM WLS? .....	23
III.1.1	Product serial number .....	23
III.1.2	DICOM C-Store Partner (SCP) .....	23
III.1.3	DICOM Worklist Partner (SCP) .....	23
III.1.4	DICOM Performed Procedure Step Partner (SCP) .....	24
III.1.5	Access rights for SIDICOM communication partners .....	24
III.2	Standard installation .....	25
III.2.1	Installing files .....	25
III.2.2	Dialog "Configure system" .....	25
III.2.3	Configuring via files .....	29
III.2.4	Creating a shortcut .....	32
III.2.5	Online help .....	32
III.2.6	Different installation directory .....	32
III.3	Client/Server mode .....	33
IV	Customizing the user interface .....	35
IV.1	Explanations concerning syntax .....	35
IV.2	Customizing the query dialog and the worklist query .....	36
IV.3	Customizing the main window .....	38

IV.4	Customizing the Performed Procedure Step Message .....	39
IV.5	Customizing the entries of the image objects (IO,DX,SC,VL,CT).....	39
IV.6	Monitoring tag values .....	40
IV.7	Setting tags to static values.....	40
IV.8	Password protection .....	40
IV.9	Overview of main configuration file "SIDICOM.ini" .....	42
V	Workflow Optimization .....	46
V.1	DICOM Worklist tag mapping .....	46
V.2	Action Item Mapping .....	46
V.2.1	What does "Action Item Mapping" mean and when is it useful? .....	46
V.2.2	Additional comments .....	47
V.2.3	Hints for intraoral images .....	47
V.2.4	Action Item Mapping configuration sheets .....	47
VI	Operational Hints (Questions & Answers, "Q & A").....	53
VII	Technical Support .....	54
VII.1	Logging of SIDICOM WLS .....	54
VII.2	Error messages.....	55
VII.3	Technical hints (Questions & Answers, "Q & A") .....	56
VII.4	Support Contact .....	58

## I Overview

### I.1 Structure of these instructions

Dear Customer,

this document is divided into four main parts:

Part I gives you a rough overview of the structure of the document, the conventions used and the main focus of the product 'SIDICOM WLS'.

Part II contains the operating instructions for users, which are intended to facilitate their daily work with 'SIDICOM WLS'.

Part III is directed at the experienced administrator or service engineer and describes installation and configuration of SIDICOM WLS. The steps for setting up the standard installation listed here have to be implemented before using the product for the first time. Here you can also find hints on how to customize configuration of the product.

Part IV lists the most common error messages of SIDICOM WLS and offers solution strategies.

## I.2 Symbols used in this manual

Key to symbols used in this manual:

Symbol	Application
<i>Italics</i>	Clicked commands in the menu bar as well as the object menus are written in italics.
"quotation marks"	The names for icons, windows, dialog boxes, pushbuttons and options are put into quotations marks.
[key]	Keys are indicated by square brackets [ ].
[key 1]+[key2]	A plus sign (+) between two keys means that both keys have to be pressed simultaneously.  Example: 'press the combination of keys [Alt]+[Esc]' means that you press the [Alt] key, hold it down and at the same time you press the [Esc] key.
	Reference to further information in another part of the manual (see contents).
	Be sure to follow this instruction.
	The following advice can facilitate working with your system and make it more simple.

### I.3 Abbreviations and Definitions

Accession number	Important hospital specific administration number in the DICOM-world.
ACR	American College of Radiology
AE	DICOM Application Entity
Client	Here: In SIDICOM Client/Server mode all those instances of SIDICOM not acting as "Server"
C-Store	DICOM-service to store objects, e.g. images
CT	Computer Tomography (DICOM image type)
DDS	Association of German dental software
DICOM	Digital Imaging and Communication in Medicine
DX	Digital X-Ray
GALAXIS	3D-SIDEXIS Image-Viewer
IO	Digital intraoral X-ray
IOD	DICOM Information Object Definition
IS	Information System (HIS, RIS, PACS within hospital)
Modalität	Image creating system (IO, Pan, ...)
Modality Worklist	See worklist
MPPS	Modality Performed Procedure Step
MWL	Modality Worklist
NEMA	National Electrical Manufacturers Association
Plug-In	Here software module to extend or change a program.
PPS	Performed Procedure Step
RawData	DICOM datatype used for raw image data or proprietary data.
SC	Secondary Capture
SCP	Service Class Provider
SCU	Service Class User
Server	Here: In SIDICOM Client/Server mode the one and only instance of SIDICOM acting as C-Store + PPS SCU
SIDEXIS	Sirona Dental X-Ray and Imaging System
SIDICOM	SIDEXIS-DICOM Interface
SLIDA	Software Link for Dental Applications
SOP	DICOM Service-Object Pair
SPS	Scheduled Procedure Step
UID	Unique Identifier, string unique in the whole network
VL	Visible Light (Photographic Image)
Worklist	Service, to get orders for the modality.

### I.4 Focus

With the help of SIDICOM WLS it is possible to integrate SDEXIS in a DICOM-compatible information and filing system. Especially the workflow can be improved considerably.

SIDICOM WLS supports:

- **Modality Worklist SCU (Service Class User)**  
Using this service enables you to retrieve order data inclusive patients' data from a DICOM compatible information system of the hospital.<sup>1</sup>
- **C-Store SCU**  
This way image data and patients' data can be sent to a DICOM system in the **SC format** (Secondary Capture), **DX format** (Digital X-Ray), **IO format** (Intraoral), **VL-format** (Visual Light), or **CT-format** (Computer Tomography). Sending DICOM Raw Data (RAW) for SDEXIS examinations is also supported.<sup>2</sup>  
It is thus possible to store the sent images in a DICOM archive or to look at them at a viewing station. This depends solely on the installed DICOM equipment.
- **Modality Performed Procedure Step (MPPS)**  
Having successfully made an exposure respectively having transferred it, a DICOM Performed Procedure Step SCP can be informed about those steps.

### I.5 System Requirements

- SDEXIS XG V.2.2 or higher is required.
- Microsoft Windows XP, recommended with Service Pack 2.
- Microsoft Windows Vista
- The installation and registration of the plug-in must be performed by an administrator or a power user.

### I.6 Scope of delivery of SIDICOM WLS

- Data carrier containing SIDICOM WLS software
- Configuration instructions/Operating instructions (this document)
- License Agreement
- DICOM Conformance Statement

---

<sup>1</sup> The respective partner must be compatible with DICOM 3.0 Modality Worklist SCP (**Service Class Provider**).

<sup>2</sup> The respective partner must be compatible with DICOM 3.0 C-Store SCP (**Service Class Provider**) and support the format in question.

### I.7 New features of the current version

In comparison to Version 2.0, the SIDICOM contains the following additional features and changes:

#### Since Version 3.0

- Switch to SDEXIS XG Plugin (better integration into SDEXIS, easier configuration without slida partners, changed order acceptance)
- New buttons for send, creating a new DICOM study (not order-dependent) and switching to the worklist
- Support of Galileos 3D volume data and views<sup>3</sup>
- Support of new DICOM data types CT, VL (visible light color images), RAW<sup>4</sup>
- Installation support by sending a demo dataset with different DICOM formats
- Tooltip function for worklist entries, with the information about the referring physician, image reason and image type
- Dose value storage in the image format DX for new Orthophos XG exposures

#### Since Version 2.2

- Possibility to send 16bit grayscale images (in addition to 8bit)
- Merging worklist data for multiple images
- Merging worklist data for multiple send partners
- Disengageable wildcard matching for "Patient's name" and "Physician's name" in the query dialog
- Internal switch of DICOM library (Therefore manual adaptation of the configuration during an update)

#### Since Version 2.1

- Can meet special networking requirements by running in Client/Server mode
- Number of query results can be limited
- Worklist item view can be sorted ascending and descending
- Additional optical representation for currently active order
- Certain dialogs can be set password-protected optionally
- Most important configuration items can be set via dialog
- DICOM tag "Software Version" can be configured
- DICOM tag "Study Description" can be configured
- Statusbar displaying last message and operating mode
- The appearance of the graphical user interface has been adapted to Sirona's (SDEXIS) neXt Generation design

---

<sup>3</sup> The respective partner must be compatible with DICOM 3.0 C-Store SCP (**S**ervice **C**lass **P**rovider) and support the CT or rather RawData format in question.

<sup>4</sup> The respective partner must be compatible with DICOM 3.0 C-Store SCP (**S**ervice **C**lass **P**rovider) and support the format in question.

- Possibility to set DICOM image tags to static values

## II Operating Instructions

"SIDICOM WLS is a separate program that enables SIDEXIS to communicate with other DICOM 3.0 partners.

SIDICOM WLS is comfortably linked to SIDEXIS via the plug-in interface.


The DICOM interface gets new orders from the worklist SCP and stores these items (Scheduled Procedure Steps) in an own "database". After activating the order, the order data is transmitted to SIDEXIS.

After pressing the send button, or if automatic sending after an exposure is activated, the image information is converted with the help of the DICOM-Interface from the SIDEXIS format to the DICOM format and sent to the DICOM Store-Partner.

During this, parts of the saved order data can be included into the image object, if the image is based on an order. Optionally, states can be transmitted to a DICOM Performed Procedure Step SCP.

### II.1 Workflow – Order based exposure

The following diagram shows the typical workflow with "SIDEXIS" and "SIDICOM WLS". If no worklist is available or if the images are made without an order, only the image exposure respectively "DICOM send" part are relevant.

#### II.1.1 Requesting order data and patient demographics

With the help of DICOM Modality Worklist Services it is possible for SIDICOM WLS to retrieve orders from the hospital's information system. These orders contain the original patient's demographics and the description of the examination to be carried through (image order).


The data transferred can then be passed on directly to SIDEXIS. A second input of the patient or order data within SIDEXIS is superfluous. On one hand you can thus avoid errors, for example spelling or transmission errors and on the other hand the workflow is faster.

### II.1.1.1 Query configuration

In certain cases it makes sense to request only a selection of order data from the information system.

By clicking the *Query/New* button using the command *Query/New*, you open a dialog box that helps you configure this selection. The configurations made there are valid for all further queries until they are changed.

The dialog may bear a  password protection. If this is the case, field contents can only be changed when the required password is entered before. Otherwise the dialog will be displayed with all fields read-protected.

**SIDEXIS DICOM WLS worklist query**

☐ Station application title:

Modality:  Physician's name:  Limit query result to:  items.

Patient's ID:  Patient's name:  Patient's name query, e.g. "Miller,J" for June and Joe Miller

Accession no.:

Procedure start date:

☐ Today

☐ Last x days:

☒ All days

☐ Free selection

Calendar: Juli 2007

	Mo	Di	Mi	Do	Fr	Sa	So
26	25	26	27	28	29	30	1
27	2	3	4	5	6	7	8
28	9	10	11	12	13	14	15
29	16	17	18	19	20	21	22
30	23	24	25	26	27	28	29
31	30	31	1	2	3	4	5

OK Cancel

It is, for instance, possible to just retrieve orders (SPSs) for a special ward or a special image type (modality). These or other limitations can result in improved handling and speed, because the system does not transfer and show orders which are not needed. The text fields should be filled out using the usual DICOM query syntax.

Please note: In order to optimize the workflow it is possible to add automatically the wildcard "\*"^\*" for name fields that are not empty and without "^". This is equivalent to searching for parts of last names with any first names.

This behaviour can be separately tuneable for "Patient's name" and "Referring Physician" in configuration file "SIDICOM.ini", see section:

```
[COMMUNICATION]
AddPhysicianAsterisk
AddPatientAsterisk
```

For the case of too many orders waiting or for too large a full response order list, the query result can be limited to n items.

It is possible to adjust parts of the "Worklist Query" dialog box to the local requirements; see Part III.

In general, a limitation of the query is not necessary, because many worklist SCPs transfer only the set of data of the respective station and the respective day. But some MWL SCPs do not recognize the limitation specified above. For those partners the query result can still be limited by setting a hard communication cut when a certain limit is exceeded, see entry "MaxItemLimitTol" in configuration file "SIDICOM.ini". Furthermore, via "SIDICOM.ini" a limit can be set for the maximum count of query result items given in the dialog. Is it exceeded by user interaction, when pressing <Ok> a warning will occur: "Sorry, the query result items' limit by configuration must not exceed ... items.". After acceptance by <Ok> the application will proceed using this limit.

Not all DICOM tags can be queried for using wildcards, thus returning multiple matching results. The DICOM standard describes where wildcard queries are allowed (F.ex. patient's name) and where not (F.ex. patient's ID).

## II.1.1.2 Repeat last query

If you click/press *Query/Repeat Last* [F6], a query is performed within the limits of the last configuration of the "Worklist Query" dialog box.

If desired, the worklist **Query/Repeat** [F6] can repeat automatically after a certain set period of time, see Part III.

For the duration of a query the user interface is blocked, symbolized by an hourglass. Depending on the speed of the systems, network link and quantity of worklist entries the waiting time can reach some minutes. If the waiting period is too long, one should limit the quantity of data which has to be transferred by executing the command *Query/New* or directly at the SCP.

## II.1.2 Selecting and activating an order

After the successful transfer of orders from the Worklist Service Class Provider to SIDICOM WLS the orders are shown in the main window.

### II.1.2.1 Selecting an order

#### II.1.2.1.1 Worklist view

The following picture shows you a worklist with a selected order.

Order	Order Time	Patient (Name, date of birth)	Patient's ID	Modality	Image Type	Accession no.
	20.07.07 12:00:00	Friedrich^Arne 29.05.79	09297264	OT		99750715
	20.07.07 12:00:00	Ballack^Michael 26.09.76	23793809	IO		02003845
	20.07.07 12:00:00	Lehmann^Jens 10.11.69	31224737	CT		00466792
	20.07.07 12:00:00	Hildebrand^Timo 05.04.79	01288080	IO		51822021
	20.07.07 12:00:00	MeierHuber^Fred 10.11.90	47550557	IO	IO	93306300
	20.07.07 12:00:00	Hans^Peter	49454944	OT		25040326
	20.07.07 12:00:00	Meier^Hans 04.03.87	50209415	CT	3D	49849553
	20.07.07 12:00:00	Lukas^Podolski 04.06.85	20678302	DX		69596292
	20.07.07 12:00:00	Kahn^Oliver 15.06.69	15006754	CT		13546521
	20.07.07 12:00:00	Klose^Miroslav 09.06.78	48113680	DX		04113786

At the beginning of each line you can see an icon which reflects the status and type of the order. Orders not yet executed are shown with the following symbols:

The icon of an X-ray tubehead symbolizes an intraoral image.

The icon of a panoramic X-ray unit symbolizes an extraoral image (Like panoramics, cephalometrics or transversal slices).

The 3D-X-ray icon symbolizes a 3D-scan.

The general X-ray icon symbolizes an image which has no further specification of image type.

An icon with an arrow ("→") underneath symbolizes that the corresponding order has not yet been activated.

To display the intraoral and extraoral icons on the screen, the image type must be known when placing the order. This image type must be mapped to the SIDEXIS image type in the configuration, ⇒Part III.

In column two you can see the time of the order, in column three the name of the patient and his/her date of birth, in column four the ID of the order, in column five the type of DICOM device, in column six the image type that is typical for the hospital or the description of the order and in column seven the accession number.

This information is only a meaningful selection of the order data intended to enable the user selecting the correct order. (⇒This displayed selection of order data can also be configured, see Part III.)

Entries in the worklist that have not been activated are automatically deleted from the list when starting a new query, except when also part of the new sets of data. Entries that are marked as finished or discontinued are removed from the list after x days. The automatic deletion after a new query as well as the number X can be configured by the sidicom.ini file ⇒Part III.

The entries of names, texts and numbers in the worklist view can be sorted. Click on the column header you want the entries to be sorted for. Sorting will take place immediately, visualized by an asterisk "\*" in front of the corresponding column's title. Clicking on this column again will invert the sort order.

### II.1.2.1.2 View Details

All data concerning a selected worklist entry transferred from the information system is listed in the lower part of the window called *View/Worklist Details*. The information is shown in the typical DICOM syntax (tag number, tag description, tag value).

The screenshot shows the SIDEXIS DICOM WLS application window. The top menu bar includes 'Query', 'Order', 'View', and 'Help'. Below the menu is a toolbar with icons for patient, order, and status. The main window is divided into two sections. The upper section displays a table of worklist entries. The lower section shows the detailed DICOM metadata for the selected entry.

Order	Order Time	Patient (Name, date of birth)	Patient's ID	Modality	Image Type	Accession no.
20.07.07 12:00:00	20.07.07 12:00:00	Friedrich^Arne 29.05.79	09297264	OT		99750715
20.07.07 12:00:00	20.07.07 12:00:00	Ballack^Michael 26.09.76	23793809	IO		02003845
20.07.07 12:00:00	20.07.07 12:00:00	Lehmann^Jens 10.11.69	31224737	CT		00466792
20.07.07 12:00:00	20.07.07 12:00:00	Hildebrand^Timo 05.04.79	01288080	IO		51822021
20.07.07 12:00:00	20.07.07 12:00:00	MeierHuber^Fred 10.11.90	47550557	IO	IO	93306300
20.07.07 12:00:00	20.07.07 12:00:00	Hans^Peter	49454944	OT		25040326
20.07.07 12:00:00	20.07.07 12:00:00	Meier^Hans 04.03.87	50209415	CT	3D	49849553
20.07.07 12:00:00	20.07.07 12:00:00	Lukas^Podolski 04.06.85	20678302	DX		69596292
20.07.07 12:00:00	20.07.07 12:00:00	Kahn^Oliver 15.06.69	15006754	CT		13546521
20.07.07 12:00:00	20.07.07 12:00:00	Klose^Miroslav 09.06.78	48113680	DX		04113786

ID: (0008,0000)	IdentifyingGroupLength =	40
ID: (0008,0050)	AccessionNumber =	02003845
ID: (0008,0080)	InstitutionName =	
ID: (0008,0090)	ReferringPhysiciansName =	
ID: (0008,1110)	ReferencedStudySequence :	
ID: (0010,0000)	PatientGroupLength =	152
ID: (0010,0010)	PatientsName =	Ballack^Michael
ID: (0010,0020)	PatientID =	23793809
ID: (0010,0030)	PatientsBirthDate =	19760926
ID: (0010,0040)	PatientsSex =	
ID: (0010,1000)	OtherPatientIDs =	
ID: (0010,1001)	OtherPatientNames =	
ID: (0010,1020)	PatientsSize =	1.89
ID: (0010,1030)	PatientsWeight =	84.0
ID: (0010,2000)	MedicalAlerts =	
ID: (0010,2110)	ContrastAllergies =	
ID: (0010,2154)	PatientsTelephoneNumbers =	
ID: (0010,2160)	EthnicGroup =	
ID: (0010,2160)	AdditionalPatientHistory =	
ID: (0010,4000)	PatientComments =	
ID: (0020,0000)	ImageGroupLength =	64
ID: (0020,0004)	StudyInstanceUID =	1.2.40.0.13.0.192.168.15.58.18481629.1147082817253.32772

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In certain cases it is possible here to view information on the patient or the order that is normally not otherwise available. Information on a pregnancy or allergies might be entered here, for example.

## II.1.2.2 Activating an order

To process an order of the worklist, it must be transferred to SIDEXIS.

The transfer can be activated through "Selecting an order" and "Order/Transmit->SIDEXIS" or pressing the "SIDEXIS" button.

Alternatively, you can transfer the order to SIDEXIS by double-clicking on it.

When the order is activated, the device icon starts blinking. Blinking stops when the next order is activated, so there is always not more than one active order blinking. In Client/Server mode there can be more than one order blinking since orders activated from other Clients are displayed as well. An activated order will only pass on to a different status (will only stop blinking), if an image concerning this order has been sent to the default partner, or if the user closes the order manually.

If the order is rejected in SIDEXIS, or if an order involves several images, it is possible to re-transfer already activated orders to SIDEXIS.

Closing an order manually might be advisable e.g., if the order should not be executed anymore, or if an order has already been executed.

In this case the two possibilities

*Order/Finish* [F8]

Order has been executed and

*Order/Discontinued* [ALT]+[F5] Order has been set discontinued

are available.

## II.1.3 Executing an order

### II.1.3.1 Accepting an order

After an order has been activated, it is transferred to SIDEXIS. If details of the patient assigned to this order do not yet exist in the SIDEXIS database, the original patient's data are automatically imported into it. A double input of patient demographics is thus avoided. If at this time an other patient is registered in "SIDEXIS", a message appears whether the patient should be checked out or the examination should be saved.

When mapping DICOM names to SIDEXIS names, the additional name designations in first and last name are adopted. This function can be switched off  
⇒ Part III "SIDICOM.ini", "MergeDICOMNames".


If not already done, the patient will be automatically registered in SIDEXIS. Within SIDEXIS the order will be accepted.

Registering the patient and getting the order automatically is only possible if - at the moment of transferring the order - SIDEXIS is in a mode where it can be interrupted. (Open dialogs as for instance during the exposure readiness are not allowed) ⇒ SIDEXIS user manual.

In the next step, SIDEXIS directly switches to image exposure (see image), if it can differ the device type on the basis of the image type and only one x-Ray component is registered. Depending on the SIDEXIS settings, the tooth number can be specified before the intraoral exposure.

If no distinction is possible, you have to select the respective type of device manually. ⇒ Chap. III)

Reasons for the order, as well as the order type, is shown in the DICOM worklist data. ⇒ Chap. III

## II.1.3.2 Image exposure with SIDEXIS (Transfer to SIDICOM WLS)

Exposures in SIDEXIS are done in the usual way. ⇒ SIDEXIS user manual.

Next, the exposure is sent via SIDICOM WLS to the desired DICOM partner (archive or viewing station). The transfer is started – depending on the configuration -

(⇒ Chap. III) directly after the exposure or after pressing the send button in SIDEXIS.

Send button in SIDEXIS.

If an automatic transfer of the exposure before image optimization is configured, it is always the **original image** (without any text on findings) that is transferred. If automatic transfer after image optimization is configured, the **optimized** image (also without any text of findings) is transferred. In case of a manual transfer one can choose between different partners. Here we transfer to the partner the **actual view of the image** together with any text on findings that may be available.

Via this send button, an image for which there is no order can be sent. In this case, however, the respective DICOM worklist data can not be integrated into the image.

SIDICOM WLS can be configured in ways that allow sending to different DICOM partners. If more than one partner exist, a dialog box appears on the screen after having pressed the send button (the default partner is preselected). There you can choose e.g. transferring the image to an archive or a viewing station. Automatic sending is only possible to the default partner.

## II.1.3.3 Sending an exposure to DICOM partner

As soon as SIDICOM WLS has received from SIDEXIS the image data to be transferred, the status of the order changes from "activated" (blinking icon) to "queued".

This is symbolized by an animated letter symbol . From now on SIDICOM WLS tries to transmit the image to the selected partner.


The transfer process takes place in the background separately from work on other orders. Depending on the configuration the user sees message boxes that indicate whether the sending procedure succeeded or failed.


Order states are only updated for the default C-Store partner. This helps in ensuring that the image is transmitted to the "main archive".

Example: notifications of whether the transfer process was successful or failed.

### Successful:

### Failed:

In the upper part of the dialog you can see the name of the patient whose image could not be transferred, the type of services (here C-Store), the partner, (here "CONQUESTSRV1"), the time of the order and the notification of successful completion or error with details of error.

These detailed dialog boxes make it possible to assign them to the respective order even if several processes of transfer are in work.


The communication or the attempt of a communication is also symbolized by an animated SIDICOM icon in the system task-bar.

### II.1.3.4 Multiple and Status Exposures

If an order is submitted that requires multiple X-ray exposures, it needs only be transmitted once to SIDEXIS. All new SIDEXIS exposures are automatically completed with the corresponding worklist data before sending them to the DICOM partner.

The order is only valid until a new order is accepted or the patient is checked out, though.

## II.2 Workflow – Creating a new DICOM study


With the "Create New DICOM Study"  command, complete SIDEXIS examinations can be sent to the archive. Here SIDEXIS creates a new DICOM study, which is not linked with a DICOM worklist order. (The images that are placed in the SIDEXIS examination can be based on different orders or created without orders.) The feedback whether a send is successful or not is given via popup dialogs.

With the default settings, the special DICOM Tags are watched during the send process (check whether information is complete) Because of that, during the send a dialog appears that gives the possibility to fill in an "Accession number". For details about tag watching, please see chapter IV.6 "Monitoring tag values"

## II.3 Link tests

To test the communication link between SIDICOM WLS and the partner, the program offers you a comfortable test function. It is called up by *Query/Echo Test...* Here, with the help of the "test" button, you can communicate directly with the DICOM partner "store", "worklist" and "PPS" via the DICOM echo service. Other non-configured partners can be directly specified with the help of the fields "Application title", "Host name" or "Port".

The test result is shown in the "status" window.

The DICOM echo service only works if the partner supports an echo SCP service. SIDICOM WLS works as a SCU.

Before you start a DICOM communication for the first time, you should ensure that the computer network is correctly configured and that the target computer can be reached via TCP/IP. Executing the "ping" command of the operating system (TCP/IP-Stack) here can be very helpful.

Note that entering Application title, Host name or Port is only possible if the echo test mode is set to "Other". In all other modes the settings are automatically retrieved from the current configuration.

## II.4 Send tests


To check which image types the archive supports, the program additionally offers a comfortable test function.

It can be called up by *Query/Send Test...*. Here it is possible to choose the different image formats and send the test-image data direct to the DICOM C-Store partner when pressing the "OK" button.

With the send test function it can be tested which DICOM-Image types the configured Store partner can receive.

Whether a link exists can be tested beforehand with the Echo-test.

Once the "OK" button is pressed, the test-image data is queued to send.

After processing the data, a positive or negative notification concerning the send process is shown.

The store partner must be configured correctly in order for the DICOM C-Store service to function properly. (Test with the echo service.)

The send Test can be used to detect which DICOM image-types in the configuration dialog can be switched on for the store-partner.

If the send process fails for certain image types, the type is not supported by the archive or it is not yet activated there.

Some archives block sending the same test objects twice. They handle a second send with the same UID as an error. In such cases, you should check if it is possible to configure the archive so that it ignores or overrides data with the same UID.

## III Installation and configuration instructions

### III.1 What is required prior to installing SIDICOM WLS?

#### III.1.1 Product serial number

The Supplement "License Agreement" provides this information.

#### III.1.2 DICOM C-Store Partner (SCP)

The following receiver information is required:

DICOM AE	DICOM <b>A</b> pplication <b>E</b> ntity Title resp. "remote_application_title" is the name of the DICOM communication partner
PortNumber	"port" is the TCP/IP port on which the remote DICOM system is listening.
HostName	"TCP/IP name" is the name of the remote host as it is known to the TCP/IP system (IP – address or alias). The address or alias may be 1 to 19 characters long (no blanks allowed!)  Example: "COMP123" or "192.168.15.7"

Furthermore the information which image formats can be processed by the DICOM C-Store SCP (image receiver). Depending on the configuration, DICOM IO (intraoral), DX (digital X-Ray), VL (visible light) and CT (computer tomography) exposure formats can be transferred. Sending RAW.Data (for SIDEXIS examinations), can be switched on or off as well. If the store partner does not support one of these formats, this data will be excluded (in the case of RAW and CT) or sent as SC (secondary capture) format instead (in the case of IO, DX and VL).

#### III.1.3 DICOM Worklist Partner (SCP)

If you want to use the worklist function, following information on the worklist partner is necessary:

DICOM AE	DICOM <b>A</b> pplication <b>E</b> ntity Title resp. "remote_application_title" is the name of the DICOM communication partner
PortNumber	"port" is the TCP/IP port on which the remote DICOM system is listening.
HostName	"TCP/IP name" is the name of the remote host as it is known to the TCP/IP system (IP – address or alias). The address or alias may be 1 to 19 characters long (no blanks allowed!)  Example: "COMP124" or "192.168.15.47"

### III.1.4 DICOM Performed Procedure Step Partner (SCP)

If you want to use the **Performed Procedure Step** function, following information on the PPS SCP is required:

DICOM AE	DICOM <b>A</b> pplication <b>E</b> ntity Title resp. "remote_application_title" is the name of the DICOM communication partner
PortNumber	"port" is the TCP/IP port on which the remote DICOM system is listening.
HostName	"TCP/IP name" is the name of the remote host as it is known to the TCP/IP system (IP – address or alias). The address or alias may be 1 to 19 characters long (no blanks allowed!)  Example: "COMP125" or "192.168.15.12"

### III.1.5 Access rights for SIDICOM communication partners

When facing communication problems, the mistake often lies in the way the access rights of the communication partners are set up. It means the partner blocks the communication with SIDICOM WLS, because the respective access rights are not set or wrongly set. Please ensure prior to the first communication that the access rights have been correctly set by the responsible administrators and that they are at hand for further questions.

The configuration of partners strongly varies; in general it is necessary to enter following information into the system:

- The **DICOM AE title** of **SIDICOM WLS**. The default entry is **SIDICOM**, but you can change it in dialog "Configure System".
- The **HOST** name or the **IP** address of the computer under which SIDICOM WLS is running. (Ask your local administrator.)
- Concerning the **C-Store** partner, it is often necessary to configure the **DICOM IODs** which are allowed to be sent. SIDICOM WLS is able to send **SC**, **VL**, **CT**, **RAW**, **IO** and **DX** IODs depending on the configuration.

## III.2 Standard installation

### III.2.1 Installing files

The recommended directory for installation is "C:\Program Files\Sirona\SIDICOM WLS", on the same computer on which images are to be sent from SIXELIS to the DICOM system via the interface.

On the installation medium, start file "Setup.exe" and follow the instructions. Thus all necessary files and subdirectories are transferred into the destination folder and the Plug-in is registered in SIXELIS.

(In the „SIXELIS Service Dialog“ the checkbox „Further communication partners“ is marked from the . This activation is necessary for an accurate order transfer and should not be deactivated manually.)

#### III.2.1.1 Update of version 2.2

The update from version 2.2 is automated as far as possible. Because the link between SIXELIS and SIDICOM WLS is changed from SLIDA to Plug-in interface, the configured DICOM partners in SIXELIS will taken over and deleted out of the SLIDA configuration. Other SLIDA partners stay unchanged.

Please install, with the help of the setup, in the (old) installation directory.

#### III.2.1.2 Updates of versions less than 2.2

An update from versions older than 2.2 is not possible. In this case, we suggest a new installation after uninstalling (deleting) the old SIDICOM version. It is recommended to save the configuration parameters first.

#### III.2.1.3 Registering the Plug-in in SIXELIS

During the installation, the SIDICOM WLS Plug-in is registered automatically.

If there are problems during the registration, the setup shows a corresponding message. In that case the Plug-in can be registered manually:

In SIXELIS choose "view/customize" , than press button "Add..." in the "Plugin Manager". Browse to the installation directory and select the file "Dipplugin.exe" (Default C:\Program Files\Sirona\SIDICOM WLS).

After the Plug-in is registered a successfully, in SIXELIS following symbols appears in the Plug-in toolbar.

The "Get DICOM Worklist" command.

The "Send to DICOM" command.


The "Create New DICOM Study" command.

### III.2.2 Dialog "Configure system"

This dialog contains the most important configuration items. Run SIDICOM WLS and activate menu command "Help/Configure system...".

The dialog may bear a  password protection. If this is the case, field contents can only be changed when the required password is entered before. Otherwise the dialog will be displayed with all fields read-protected.

**Configure system**

**Common**

DICOM Application Entity Title (AE Title) SIDICOM

Institution name No Institution Name

Station name No Station Name

Product Serial No. 0

SIDICOM data path (Client/Server mode) >>

**Mode**

☒ Stand-alone

☐ Client (No direct C-Store SCP contact)

☐ Server (Direct C-Store SCP contact only here)

**Storage SCP**

DICOM Application Entity Title (AE Title)

Host name

Port number 0

☒ Supports DICOM ID format

☒ Supports DICOM DX format

☒ Supports DICOM VL format

☒ Supports DICOM CT format

☒ Supports DICOM RAW format

☐ Automatically receives a copy of new exposures

☐ Before image optimization

☒ After image optimization

**Modality Worklist SCP**

DICOM Application Entity Title (AE Title)

Host name

Port number 0

☐ Query automatically every [seconds] 60

Insert further SCP's by editing the corresponding .ini files manually (Refer to manual).

OK Cancel

The data path for the client server mode can be browsed for by pressing the adjacent ">>" button.

**In "Stand-alone" mode (default) the "SIDICOM data-path" is empty!**

Just leave all fields of items you do not need blank.

Following some comments to the fields are listed:

Field	Comment
<b>Common</b>	<b>Common items for SIDICOM WLS</b>
DICOM Application Entity Title (AE Title)	DICOM Application Entity title of SIDICOM WLS (Normally "SIDICOM")
Institution name	Insert here the name of your institution
Station name	Insert here the name of the ward
Product Serial No.	Sirona serial number of SIDICOM WLS. The serial number can be found in the supplement of the "License Agreement", as well. It is issued by Sirona. (This number has only numerical digits.)
Data path	In Client/Server mode, insert the path where this instance will place, and look for, respectively, commonly used data, like storage requests. This directory should be common to every client and to the server. Clients as well as the server require read and write access to this directory.  In "Stand-alone" mode the "SIDICOM data-path" is normally empty.
<b>Mode</b>	<b>The mode in which this instance of SIDICOM WLS shall run</b>
Stand-alone	Standard, no Client/Server mode (Default)
Client (No direct C-Store SCP contact)	Client/Server mode, act as a Client
Server (Direct C-Store SCP contact only here)	Client/Server mode, act as the (one and only!) Server
<b>Storage SCP</b>	<b>DICOM settings for the (first/default) C-Store SCP</b>
DICOM Application Entity Title (AE Title)	Insert here this partner's DICOM Application Entity title
Host name	Insert this partner's TCP/IP name here
Port number	Insert this partner's TCP/IP port number her
Supports DICOM IO format	If checked, SIDICOM WLS assumes that this C-Store SCP accepts (and that viewers in use elsewhere can read) DICOM IOD "IO". Otherwise, SIDICOM WLS will send such images as "SC" (Secondary Capture) objects instead.
Supports DICOM DX format	If checked, SIDICOM WLS assumes that this C-Store SCP accepts (and that viewers in use elsewhere can read) DICOM IOD "DX". Otherwise, SIDICOM WLS will send such images as "SC" (Secondary Capture) objects instead.
Supports DICOM VL format	If checked, SIDICOM WLS assumes that this C-Store SCP accepts (and that viewers in use elsewhere can read) DICOM IOD "VL". Otherwise, SIDICOM WLS will send such images as "SC" (Secondary Capture) objects instead.
Supports	If checked, SIDICOM WLS assumes that this C-Store SCP accepts (and that

Field	Comment
DICOM CT format	viewers in use elsewhere can read) DICOM IOD "CT". Otherwise, SIDICOM WLS will <b>not send</b> such images.
Supports DICOM RAW format	If checked, SIDICOM WLS assumes that this C-Store SCP accepts (and that viewers in use elsewhere can read) DICOM IOD "RAW". Otherwise, SIDICOM WLS will <b>not send</b> such images.
Automatically receives a copy of new exposures	This means that the C-Store-SCP receives the images directly after exposure. If checked, it must be defined whether the copy is send before or after image optimisation. This setting has only effect for the configured main archive under "Storage SCP".
<b>Modality Worklist SCP</b>	<b>DICOM settings for the Modality Worklist SCP</b>
DICOM Application Entity Title (AE Title)	Insert here this partner's DICOM Application Entity title
Host name	Insert this partner's TCP/IP name here
Port number	Insert this partner's TCP/IP port number here
Query automatically every	Automatic worklist query on?
[seconds]	Duration of the automatic worklist query cycle [s]

As in the stand-alone case and for the Server in Client/Server mode, you can leave the "Data path" field blank. SIDICOM WLS will then automatically choose its installation directory as "Data path". Of course, in Client/Server mode, the Clients still need access to this path. A quite elegant way of specifying "Data Path" is it to enter the database of SIDEXIS here, which must be already known to all workstations with SIDEXIS installed.

In later configuration sessions, be sure not to select the wrong "Data path". SIDICOM WLS will use the directory given here to create a subdirectory structure of its own there, normally starting with "\SIDICOM\". So, the field "Data path" should never end with "\SIDICOM\". Otherwise, SIDICOM WLS may loose data already stored in the right place.

### III.2.3 Configuring via files

It may be necessary to adapt some files with the help of a pure text editor, e.g. Notepad. This procedure is described in the following

#### III.2.3.1 SIDICOM configuration file

##### SIDICOM.ini

Default in directory "C:\Program Files\Sirona\SIDICOM WLS" or "C:\SiDICOM" in the case of update.

##### **[Hostname]**

ProcedureStep=     ⇒ Insert the DICOM application entity title of the MPPS SCP here

*If you do not have a worklist or PPS partner, leave a blank field (as shown in the example of ProcedureStep). If you work with several C-Store partners a choice box appears when sending images manually.*

Every default configuration partner (Store, Worklist or PPS) runs its own section for the special association parameters. In the case of a Worklist or PPS partner the **section** is named as the AE-title of the partner. **In the case of a C-Store partner since V.3.0 the name of the partner section is build with the AE-title of the C-Store partner and a "CSTORE\_" prefix.**

Section sample for a C-Store partner with the AE-Title "CONQUESTSRV1":

[CSTORE\_CONQUESTSRV1]

Section sample for a PPS partner with the AE-Title "PPSSERV":

[PPSSERV]

The sections for the Worklist partner and for the default Store partner can be generated and changed via dialog "Configure system".

You should create a separate section for each additional C-Store partner (which can be an archive as well as a viewing workstation). In addition to the association parameters you can set up here which kind of DICOM IODs (Image types) are supported by the partner. If "SupportIO", "SupportVL" or "SupportDX" is equal to "0", the images are sent as "SC" (Secondary Capture) objects instead of the "IO" (Intraoral), "VL" (Visible light) and "DX" (Digital X-Ray) image types.

In the case "SupportCT" or "SupportRAW" equal to "0", no "RAW" (DICOM RAW) or "CT" (Computer Tomography) data is transferred.

##### **C-Store partner example**

Generate for each store partner a separate section with its own AE-Title as name.

[CSTORE\_ARCHIVE2]

SupportIO=1

SupportDX=1

SupportVL=1

SupportCT=1

SupportRAW=1

##### **Example C-Store AE**

⇒ Set the entry to "1", if the partner supports the IO format

⇒ Set the entry to "1", if the partner supports the DX format

⇒ Set the entry to "1" if the partner supports the VL format

⇒ Set the entry to "1" if the partner supports the CT format

⇒ Set the entry to "1" if the partner supports the RAW format

HostName=192.168.15.01

⇒ Hostname or IP-address of the partner

PortNumber=59550

⇒ Port number of the partner

Profile=StoreProfileDefault

⇒ Link to a communication profile in "Profiles.cfg"

*The Mailslot files must be the same files as of the respective partner described in the chapter "Configuration of SIDEXIS".*

### Worklist partner example:

[WORKLIST1]

Example Worklist AE

HostName=192.168.15.02

⇒ Hostname or IP-address of the partner

PortNumber=59551

⇒ Port number of the partner

Profile=WorkPPSPProfile

⇒ Link to a communication profile in "Profiles.cfg"

### PPS partner example:

[PPS1]

Example PPS AE

HostName=192.168.15.03

⇒ Hostname or IP-address of the partner

PortNumber=59552

⇒ Port number of the partner

Profile=PpsProfileDefault

⇒ Link to a communication profile in "Profiles.cfg"

### [SLIDA]<sup>5</sup>

In the "SLIDA" section you can change the general settings for the DICOM-Tag transformation from SIDICOM to SIDEXIS. The presets are:

PermanentDentist=<(0032,1032)> ;Requested Physician

ImageReason=<(0040,0100)(0040,0007)> ;Scheduled Procedure Step Description

ImageType=<(0040,0100)(0040,0008)(0008,0100)> ;Scheduled Action Item Code Value

"Requesting physician" will be mapped to the "Permanent Dentist" in SIDEXIS. "ImageReason" and "ImageType" are the names that SIDEXIS is using for the orders. They can be mapped from any DICOM tag.

If you move the mouse slowly over a worklist item, a tool tip appears. Here the contents of the fields "PermanentDentist", "ImageReason" und "ImageType" are shown. In the example, only the fields "ImageType" and "ImageReason" are engaged.

---

<sup>5</sup> Even if SIDICOM V.3.0 or higher is not longer based on "Slida", these entries are still used to be compatible with older versions.

SIDEXIS DICOM WLS							
Query Order View Help							
Order	Order Time	Patient (Name, date of birth)	Patient's ID	Modality	Image Type	Accession no.	
	08.05.06 13:07:50	Friedrich^Arne 29.05.79	09297264	IO	IO	99750715	
	08.05.06 13:11:40	Ballack^Michael 26.09.76	23793809	IO	IO	02003845	
	08.05.06 11:12:48	Lehmann^Jens 10.11.69	31224737	PX	PX	00466792	
	08.05.06 10:27:33	Hildebrand^Timo 05.04.79	01288080	CT	CT	51822021	
	04.10.05 10:25:00	MeierHuber^Fred 10.11.90	47550557	IO	IO	93306300	
	04.10.05 10:24:58	MeierHuber^Fred 10.11.90	47550557	DX	DX	49849553	
	04.10.05 10:24:58	MeierHuber^Fred (RequestingPhysician)	593	CT	CT	58559456	
	08.05.06 11:06:58	IO 1,1 Zahnaufnahme (ScheduledProcedureStepDescription)	13546521	DX	DX	13546521	
	04.10.05 10:25:02	Baumgart^Inge 05.09.78	63877561	IO	IO	10907933	
	05.09.07 13:01:55	MeierHuber^Fred 10.11.90	47550557	CT	CT	93306300	
	08.05.06 13:21:31	Klose^Miroslav 09.06.78	48113680	PX	PX	04113786	
	08.05.06 13:14:15	Podolski^Lukas 04.06.85	20678302	DX	DX	69596292	
	04.10.05 10:25:01	Schweitzer^Malis 02.03.55	88465547	OT	OT	38474277	
ID: (0008,0000) IdentifyingGroupLength = 40 ID: (0008,0050) AccessionNumber = 93306300 ID: (0008,0080) InstitutionName = ID: (0008,0090) ReferringPhysiciansName = ID: (0008,1110) ReferencedStudySequence : ID: (0010,0000) PatientGroupLength = 152 ID: (0010,0010) PatientsName = MeierHuber^Fred ID: (0010,0020) PatientID = 47550557 ID: (0010,0030) PatientsBirthDate = 19901110 ID: (0010,0040) PatientsSex = ID: (0010,1000) OtherPatientIDs = ID: (0010,1001) OtherPatientNames = ID: (0010,1020) PatientsSize = ID: (0010,1030) PatientsWeight = ID: (0010,2000) MedicalAlerts = ID: (0010,2110) ContrastAllergies = ID: (0010,2154) PatientsTelephoneNumbers = ID: (0010,2160) EthnicGroup = ID: (0010,2166) AdditionalPatientHistory =							
* 13 worklist items found based on the criteria entered. Ident: 'MKC01' Info: 'Worklist Query' OrderTime: '16:27:09'							Stand-alone

Since many hospitals often have their own names for image types field "Image Type" is not derived directly from a DICOM tag but is rather mapped via a relational table ⇒ [ActionItemMap].

In case you already know the image type when placing the order you can map the institute's internal name for the image type via an entry in the "ActionItemMap" to get a SIDEXIS-compatible image type.

Please read chapter "Action Item Mapping" for further details!

The mapping syntax is described in Chapter IV.1.

In order to deploy the new DICOM image format "IO" it is necessary to store the area of examination or the tooth number within the image. Ideally the position is already known when placing the order or it is defined by SIDEXIS before or after the exposure. Please note that when sending the image directly after the exposure automatically, there is no further possibility to set the image type.

### III.2.3.2 DICOMcommunication profile ("Profiles.cfg")

In file "Profiles.cfg" in the SIDICOM installation directory, there are preconfigured communication profiles of the partners. Normally the file need not be adjusted. Please change this file only upon advice of the Sirona Customer Service or if you are familiar with DICOM communication protocols.

Exception: In order to communicate with the DICOM server from Schick Technologies Inc. it is necessary to reduce the transfer syntaxes to one, exactly "LittleEndianImplicit". To do that please change section "TransferSyntax\_Uncompressed" as follows.

[TransferSyntax\_Uncompressed]  
TransferSyntax1=LittleEndianImplicit

Remove all additional transfer syntaxes ("TransferSyntax2",3,...).

### III.2.4 Creating a shortcut

For an additional starting possibility of SIDICOM WLS you can create a shortcut on the desktop or generate an entry in the Windows system menu.

To achieve this, there are a lot of possibilities, which have effects on the actual user or all users. Please look in the Windows help for further details.

Setting up a shortcut for the actual user on the desktop:

Open the Explorer.

Select the file to be connected ("C:\Program Files\Sirona\SIDICOM WLS\Sidicom.exe").

Now click the right mouse key.

Select the entry "Copy" in the context menu.

Now place the mouse pointer on the desktop and press the right mouse key.

Select the entry "Create shortcut" in the context menu.

### III.2.5 Online help

The English version of this document is also used as the online help of SIDICOM WLS.

If you press "F1" or *Help/Help topics*, the document "Manual.pdf" in the SIDICOM installation directory is shown. To read the document, it is necessary that the Acrobat Reader® is already installed. This software is freely available at <http://www.adobe.com/acrobat/>. The English and German version can be found on the shipped CD too.

If you want to use the German version of this document as the online help, please replace file "Manual.pdf" with the German one. (This shall be renamed in "Manual.pdf".)

### III.2.6 Different installation directory

To choose a different installation path for a new installation, please proceed as follows.

On the installation medium, start the file "Setup.exe" and follow the instructions. Thus all necessary files and subdirectories are transferred into the destination folder.


Please assure that the users have the required rights (full access) for the installation directory.

## III.3 Client/Server mode

This mode is to be seen optional as an alternative to the normal Stand-alone mode. This is necessary when, for example due to license requirements, only one workstation shall use the C-Store or MPPS service.

Following scheme visualizes a possible scenario. The different symbols can, but need not, be implemented with dedicated hardware. One of the workstations where SIDICOM WLS is running, at the same time is the Server for the other two, i.e. it takes over for itself as well as for Client A + B the complete communication with the C-Store and MPPS SCP. By contrast every SIDICOM WLS can query the Worklist SCP. The Clients are done processing an order when they have deposited the files created and to be transferred on the File Server, where the "Data path" was published. There the Server takes them up in the set intervals for further completion. However the Clients still can display the current degree of completion, since they have access on that via "Data path", too. In the intervals set for them they update the list displayed locally at their site which then always accounts for the whole state of all orders in the network.

Changes in the worklist data (SPSs), for example through a new query or through a state change at one instance, thus will affect the displays of ALL Clients AND of the Server.

The Refresh of the display takes place during idle time within a configurable time interval, otherwise before and after every important action. For example, an image is sent by the Server; the Clients realize that at the latest after the interval, or Client A does a new query, Client B takes over the result of the new query at the latest after the interval as well automatically. On this occasion resorting will take place if sorting was in effect before – items will always be unsorted after queries.

Reports concerning success or failure of data sent to C-Store or PPS SCP are displayed only at the Server, though.

A further "though": Local workstation timestamps may differ from each other. Storage requests are handled according to their relative timestamp.

In order to optimize quick Server responses to store requests by the Clients, assure by some decent, standard means that the clocks of all workstations running SIDICOM instances, be it Client or Server, be synchronized.

## IV Customizing the user interface

To make SIDICOM WLS more convenient for our customers, it is possible to configure file "SIDICOM.ini" and adapt many settings to local needs. "SIDICOM.ini" partially represents the configuration interface of SIDICOM WLS. For example different DICOM tags can be mapped to SIDEXIS data or can be shown in SIDICOM WLS.

### IV.1 Explanations concerning syntax

- Text without brackets is a normal text.
- A comment begins with a semicolon ";".
- A DICOM tag is written in parentheses as in standard use. The group and the element are separated by a comma, both having four characters always. Example: (0010,21C0) is the tag for Pregnancy Status
- Pointed brackets "<>" enclose a variable. It could be a DICOM tag or a SIDICOM variable. They are used in assignments, when you have to differentiate between a direct string as value and the contents of a tag, or when you have to indicate a certain field within the surface of SIDICOM WLS.
- To identify a tag of one sequence, the complete sequence path has to precede the tag. Example:
  - [SLIDA]  
ImageType=<(0040,0100)(0040,0008)(0008,0100)> ;Scheduled Action Item Code Value
  - Here the SIDEXIS/Slida variable ImageType is described with the value of tag Scheduled Action Item Code Value (0008,0100). Since this tag is embedded in the "Scheduled Action Item Code Sequence" (0040,0008) and this sequence is embedded in the "Scheduled Procedure Step Sequence" (0040,0100) of a worklist message, the whole sequence path is indicated.
  - If you generate your own messages, use following syntax:
  - Example: Extract from the worklist query message.
  - [WQuery]  
(0010,0030);Patient's Birth Date  
<Sequence>  
REF\_STUDY  
(0008,1110);Referenced Study Sequence  
(0008,1150);Referenced SOP Class UID  
(0008,1155);Referenced SOP Instance UID  
<!Sequence>

Here the worklist query message is composed. It contains the tag (0010,0030) and a sequence with the name REF\_STUDY. The sequence contains the tags (0008,1110), (0008,1150) and (0008,1155). <Sequence> identifies the beginning of the sequence, <!Sequence> the end. Sequences can be embedded in other sequences.

## IV.2 Customizing the query dialog and the worklist query

Only change these settings, if you know the DICOM standard well. Wrong entries can lead to malfunction or data errors. If in doubt, please contact the Sirona Hotline.

The query dialog of SIDICOM WLS offers you many ways of entering data. Sometimes you might like to use other entry fields. For this reason, you can freely configure the five text fields of the dialog (see picture).

**The headings of the edit fields** are entered in section "QueryDlg". The contents of the fields are always treated as abstract variables 2-6. Variables 1 and 8 are the station application title and the procedure start date and cannot be changed.

[QueryDlg]

UserDefine1=Modality:

UserDefine2=Physician's Name:

UserDefine3=Patient's ID:

UserDefine4=Patient's Name:

UserDefine5=Accession Number:

When issuing a DICOM worklist request, according the DICOM standard only those tags or sets of data are transferred that are requested by the SCU from the SCP. In file "SiDICOM.ini" this worklist request can be customized in detail, enabling the program to extend the data to be transferred or to limit it. In most of the cases the default can remain unchanged.

**The contents of the edit fields** are filled into different tags of the Worklist Find message and can be configured via section "WLQuery" in "SIDICOM.ini" as described in extracts below:

[WLQuery]

(0010,0010)<5>;Patient Name

(0010,0020)<4>;Patient ID

(0010,1000);Other Patient's ID

<Sequence>

```
SCHEDULED_PROCEDURE_STEP  
(0040,0100);Scheduled Procedure Step Sequence  
(0040,0001)<1>;Scheduled Station AE Title  
(0008,0060)<2>;Modality  
<!Sequence>  
...
```


Line 1 determines that only the sets of data of certain patients are to be transferred. Here the <5> marks the internal dialog variable of the query dialog pos.5. **The numbering counts here inclusive the non-configurable elements of the dialog; the patient's name is therefore pos. 5, not 4.** (If you have changed the element UserDefine4 in section [QueryDlg], you are not allowed to use its variable (var. 5) for the tag (0010,0010) Patient Name, but only for the right tag).

In line 2 the patient ID string is transferred. With line 3 you request the Other Patient's ID tag. Then the scheduled procedure step sequence, which has also to be transferred, is defined. Integrated in the sequence is e.g. the tag modality, which has the dialog variable <2>.

In any case, the DICOM type 1 tag should be requested and transferred by the SCP. The Requested Procedure ID, Scheduled Procedure Step ID and Study Instance UID are necessary for a definite reference of the SPS (Scheduled Procedure Steps).

## IV.3 Customizing the main window

Order	Order Time	Patient (Name, date of birth)	Patient's ID	Modality	Image Type	Accession no.
	08.05.06 13:07:50	Friedrich^Arne 29.05.79	09297264	IO		99750715
	08.05.06 13:11:40	Ballack^Michael 26.09.76	23793809	IO		02003845
	08.05.06 11:12:48	Lehmann^Jens 10.11.69	31224737	PX		00466792
	08.05.06 10:27:33	Hildebrand^Timo 05.04.79	01288080	IO		51822021
	04.10.05 10:25:00	MeierHuber^Fred 10.11.90	47550557	DX	IO	93306300
	12.07.07 10:57:38	Hans^Peter	49454944	OT		25040326
	04.10.05 10:24:58	Meier^Hans 04.03.87	50209415	PAN	PAN	49849553
	08.05.06 13:14:15	Lukas^Podolski 04.06.85	20678302	DX		69596292
	08.05.06 11:06:57	Kahn^Oliver 15.06.69	15006754	PAN		13546521
	08.05.06 13:21:31	Klose^Miroslav 09.06.78	48113680	PX		04113786

\*10 worklist items found based on the criteria entered.' Ident: 'RISSERV' Info: 'Worklist Query' OrderTime: '09:34:49'

Stand-alone

The column entries of the main window can be retrieved from any of the tags of a worklist entry. (Only the order column can not be changed, also the number of columns is fixed.)

Corresponding to the entries in section [WIView] here e.g. the heading "modality" and the contents of the modality tag of the Scheduled Procedure Step Sequence are mapped to the fourth column "<4>". The heading "accession number" and the contents of the respective tag are mapped to the sixth column "<6>".

The key words "string" and "number" have an influence on how the entries are sorted. (Sorting according to dates is not supported.) As you can infer from the entries of the second column "<2>", it is also possible to assign the contents of two tags to a column. Here, it is (0010,0010) for the patient's name and (0010,0030) for the date of birth. (As separators between several tags, the system always inserts a blank.)

### [WIView]

```
<1>"Order Time"<(0040,0100)(0040,0002)(0040,0100)(0040,0003)>String
<2>"Patient (Name, date of birth)"<(0010,0010)(0010,0030)>String
<3>"Patient's ID"<(0010,0020)>Number
<4>"Modality"<(0040,0100)(0008,0060)>String
<5>"Image Type"<(0040,0100)(0040,0008)(0008,0100)>String
<6>"Accession no."<(0008,0050)>Number
```

## IV.4 Customizing the Performed Procedure Step Message

Change these settings only, if you know well the DICOM standard. Wrong entries might lead to malfunctions or distortion of data. If in doubt, please contact the Sirona Hotline.

SIDICOM WLS is capable of informing a Performed Procedure Step SCP on the executed activities. Some of the tags of the PPS-Create Message are determined by worklist entries and can be influenced via the section PPSCreate of the configuration file.

The entries have a similar structure as the entries of the [WIQuery] section described above. Except that here some of the variables of the worklist are mapped to those of the PPSCreate. Excerpt from the PPSCreate section:

```
[PPSCreate]
(0010,0010)<(0010,0010)>;Patient Name
(0010,0020)<(0010,0020)>;Patient's ID
(0010,0030)<(0010,0030)>;Patient's Birth Date
(0010,0040)<(0010,0040)>;Patient's Sex
(0008,1120);Referenced Patient Sequence
...
```

E.g. item "Patient Name" of the worklist item here is mapped to the same tag of the PPS item. The same applies for the patient's ID, birthday and patient's sex.

## IV.5 Customizing the entries of the image objects (IO,DX,SC,VL,CT)

Change these settings only, if you know well the DICOM standard. Wrong entries might lead to malfunctions or distortion of data. If in doubt, please contact the Sirona Hotline.

Like the PPS, tags of the image IOD can also be allotted with information from the worklist entries or with fix values.

Extract from the ImageHeader section:

```
[ImageHeader]
(0010,0010)<(0010,0010)>;Patient Name
(0010,0020)<(0010,0020)>;Patient ID
(0010,0030)<(0010,0030)>;Patient's Birth Date
(0010,0040)<(0010,0040)>;Patient's Sex
<Sequence>
REF_PATIENT
(0008,1120)<(0008,1120)>;Referenced Patient Sequence
(0008,1150)<(0008,1120)(0008,1150)>;Referenced SOP Class UID
(0008,1155)<(0008,1120)(0008,1155)>;Referenced SOP Instance UID
<!Sequence>
```

## IV.6 Monitoring tag values

Within some DICOM environments it is observed that single tags are not or only under certain circumstances available, but these data are absolutely necessary for a later reference in the archive.


In this case you can use the monitoring function of SIDICOM WLS. This function can also be interesting, if there is no connection to a worklist.

For monitoring tags, you can enter the tags to be monitored in the section "CheckImageHeader" in the sidicom.ini file. In the following example the date of the patient's birthday and the accession number are monitored.

```
[CheckImageHeader]
(0010,0030);Patient's Birth Date
(0008,0050);Accession Number
```

Note that these tag monitor values are also the factory setting!

Prior to sending the image, SIDICOM WLS checks whether the tags to be monitored are filled in or not. If empty the following dialog for entering the missed tag value pops up:

## IV.7 Setting tags to static values

By analogy to tag monitoring section [CheckImageHeader] can be used to set static values in DICOM image tags.

For this the value is directly written after the tag.

Example:

```
(0008,0060)CT;Modality
```


Note that the value must be in the range granted to it in the DICOM standard. It is not possible to change the Modality for CT imagetypes and RawData objects. (Whether or not the value is inside the allowed range according to the standard is not tested by the program.)

## IV.8 Password protection

As an option, certain dialogs can be password-protected. When the protection is switched on, the user is asked to enter the password before the dialog wanted appears: This can be advantageous for sites where, for example, queries are configured once for all specific for the modality.

Pressing <OK> without a password entered will result in displaying the next dialog with all fields read-protected.

Ask the Sirona hotline for how to set up the password protection.

With the correct password entered once password protection is disabled for the whole rest of the current SIDICOM WLS session. In order to re-enable protection just restart the application.

## IV.9 Overview of main configuration file "SIDICOM.ini"

The following table may contain only those entries that can become subject to manual editing by the user or administrator.

Section (in brackets "[ ]") or entry	Description	Entries at delivery
<b>[Hostname]</b>	<b>Section for the DICOM AE title of every partner</b>	
Store	Name of the standard partner for the C-Store service. For this partner the worklist states are updated.	""
Worklist	Like Store the name of the worklist partners (SCP) If you do not wish a link to a worklist partner, leave a blank field.	""
ProcedureStep	Like Store the name of the Performed Procedure Step partners (SCP) If you do not wish a connection to a PPS, leave a blank field.	""
<b>[CSTORE_TIANI_ARCHIV E]</b>	<b>Sample settings for a C-Store partner (Here: Server "TIANI_ARCHIVE")</b>	
SupportIO	The partner supports the DICOM IO image format. If you enter a "0" here, the SC images are transferred.	
SupportDX	The partner supports the DICOM DX image format. If you enter a "0" here, SC images are transferred.	
SupportVL	The partner supports the DICOM VL image format. If you enter a "0" here, SC images are transferred.	
SupportRAW	The partner supports the DICOM RAW data format. If you enter a "0" here, no RAW-data objects are transferred.	
SupportCT	The partner supports the DICOM CT image format. If you enter a "0" here, no 3D images are transferred.	
HostName	Host name or IP address for the partner	
PortNumber	Port where the partner offers the service	
Profile	Communication profile of the partner (reasonably "StoreProfileDefault").	
<b>[TIANI_WORKLIST]</b>	<b>Sample settings for a Worklist SCP (Here: Server "TIANI_WORKLIST")</b>	
HostName	Host name or IP address for the partner	
PortNumber	Port where the partner offers the service	
Profile	Communication profile of the partner (reasonably "WorkProfileDefault")	
<b>[TIANI_PPS]</b>	<b>Sample settings for a PPS SCP (Here: Server "TIANI_PPS")</b>	
HostName	Host name or IP address for the partner	
PortNumber	Port where the partner offers the service	
Profile	Communication profile of the partner (reasonably "PpsProfileDefault")	
<b>[SLIDA]</b>	<b>Settings for the "SLIDA" interface</b>	

Section (in brackets "[ ]") or entry	Description	Entries at delivery
PermanentDentist	SIDEXIS permanent dentist: can be mapped from any worklist tag.	"<(0032,1032)> ;Requested Physician"
ImageReason	SIDEXIS image reason: can be mapped from any worklist tag.	"<(0040,0100)(0040,0007)> ;Scheduled Procedure Step Description"
ImageType	SIDEXIS image type: can be mapped from any worklist tag. Anyhow this tag is not taken directly from the modality tag, but is being converted according the [ActionItemMap]	"<(0040,0100)(0040,0008)(0008,0100)> ;Scheduled Action Item Code Value"
SlidInPolling	Starts an automatic polling of the SlidIn Mailslots. This is necessary, if the images have to be transferred directly after the exposure is done.	1
SlidInPollingTime	Duration of the automatic polling [s].	10
MergeDICOMNames	<p>If the default is set to "1", any name prefixes, -suffixes and middle names are integrated in the SIDEXIS first and last names before the transfer to SIDEXIS.</p> <p>An example from the Standard: Rev. John Robert Quincy Adams, B.A. M.Div. "Adams^John Robert Quincy^^Rev.^B.A. M.Div." [One family name; three given names; no middle name; one prefix; two suffixes.]</p> <p>becomes "Adams Rev. B.A. M.Div.,John Robert Quincy" [last name],[first name]</p> <p>If the entry set is to "0", only first and last names are converted.</p>	1
TransmitPatID	<p>=1 =&gt; Patient ID is transmitted to SIDEXIS</p> <p>=0 =&gt; ID is not transmitted (Useful when the Ids in the MWL may come from different management systems whose ID number ranges could possibly overlap)</p>	1
<b>[COMMON]</b>	<b>General information about SIDICOM</b>	
ApplicationTitle	DICOM AE title of the interface	"SIDICOM"
InstitutionName	Name of the institute	"No Institution Name"
StationName	Name of the ward	"No Station Name"
CharSet	Used character set	"ISO_IR 100"
SerNo	Serial number of SIDICOM	0
Validate	Conformance check of DICOM format	0
ErrorTime	Length of time for which an error message appears. 0 means that it remains until acknowledgement	0
MessageTime	Length of time for which an error message appears. 0 means that it remains until acknowledgement	20
ShowStoreMsg	Showing all messages of the Store Services?	1

Section (in brackets "[ ]") or entry	Description	Entries at delivery
ShowWorklistMsg	Showing all messages of the Worklist Service? If set to "0" no Worklist messages are displayed.	1
ShowPpsMsg	Showing all messages of the PPS Service?	0
Resolution	Do you want to transfer the resolution information tag 0028 0030 in SC image?	1
XRayPhotoInter	DICOM Photometric Interpretation. (When "MONOCHROME1" the image is being inverted.)	"MONOCHROME2"
ModalityKind	Default value of the modality type of the SC-images. In general, this value is overwritten with the modality type of the worklist.	"OT"
StudyDescription	Default value of the study description. In general, this value is overwritten with the value of the worklist.	"SIDEXIS study = images of identical date"
SWVersion	Default value for DICOM tag "Software Version", being transferred when present alternatively to the corresponding SIDEXIS field	""
WorklistQuery	Automatic worklist query on?	0
WorklistQueryTime	Duration of the automatic worklist query cycle [s]	60
ResponseTimeOut	Timeout-duration of a DICOM response.	30
SetDeliveryMax	Max number of repeated delivery attempts when having communication problems via the network.	3
SetDeliveryPause	Time between two repeated delivery attempts. [s]	90
AutoTransmit	Enable/Disable automatic transfer of a new order to SIDEXIS (The order is only transmitted if exactly one (1) new order was transferred by the worklist query)	0
GarbageDays	Duration after which number of days completed or discontinued orders are deleted from the worklist view. It is valid to use floating point values; E.g. "0.01" is equivalent to 14.4 minutes.  Recommended to set to "0.5" in most cases!	14 0.5 beginning with a new installation of version 2.2
DelOldQueryResults	Delete old, not yet started orders when doing new worklist query?	1
DataRefreshTime	Interval [s] after which the local worklist view is being updated with the centralized data – relevant in Client/Server mode only	60
ExeMode	Execution mode: 0 = Stand-alone 1 = Client 2 = Server	0
pw	Password for opening of otherwise thus protected dialogs; active when entry exists, entry in clear	""
<b>[COMMUNICATION]</b>	<b>Communication parameters</b>	
AddPatientAsterisk	=1 => Add **^** in "Patient's name" for a worklist query	1
AddPhysicianAsterisk	=1 => Add **^** in "Ref. Physician" for a worklist query	0
<b>[DEBUG]</b>	<b>Additional setting for DEBUG purposes</b>	
Debug	Special Debug mode. Additional edits and error messages to come closer to problems?	0

Section (in brackets "[ ]") or entry	Description	Entries at delivery
Verbose	Additional logging (particularly internal states)	0
LogFile	Name of logfile (path of the logfile). The default when having no entry is "log\sidicom.log" file	""
LogFileLen	Maximal length of the logfile [kByte]	1000 5000 beginning with a new installation of version 3.0
<b>[QueryDlg]</b>	<b>Configuration of the Query Dialog see III.4.3</b>	
UserDefine1	Text of the first "UserDefine" field	"Modality:"
UserDefine2	Text of the second "UserDefine" field	"Physician's Name:"
...		
MaxItemUpperBound	Maximum SPS count allowed to get during the query (Limitation for corresponding dialog item)	80
MaxItemLimitTol	Count of additional received SPS after "MaxItemUpperBound" has already been reached, when exceeded cut the connection to the MWL SCP	20
<b>[WIQuery]</b>	<b>Section to define the worklist query request message.</b>	
Special syntax	See "Customizing the query dialog and the worklist query"	Completely pre-set by factory
<b>[WIView]</b>	<b>Section to configure the worklist view window. (main window)</b>	
Special syntax	See "Customizing the main window"	Completely pre-set by factory
<b>[PPSCreate]</b>	<b>Section to define the PPS create message.</b>	
Special syntax	See "Customizing the Performed Procedure Step Message"	Completely pre-set by factory
<b>[CheckImageHeader]</b>	<b>Section to define the check of the image header tag.</b>	
Special syntax	See "Monitoring tag values" and "Setting tags to static values"  Example: [CheckImageHeader] (0008,0060)MKG;Modality In the example Tag "Modality" contains "MKG".	Completely pre-set by factory
<b>[ImageHeader]</b>	<b>Section to take map worklist data in the image header.</b>	
Special syntax	See "Customizing the entries of the image objects (IO,DX,SC,VL,CT)"	Completely pre-set by factory
<b>[ActionItemMap]</b>	<b>Mapping table to convert clinic specific image types into SIDEXIS image types</b>	
Special syntax	See "Action Item Mapping"	<b>Must be individually adjusted!</b>
<b>[POSITIONS]</b>	<b>Window &amp; column positions of the SIDICOM WLS's main window</b>	
div. entries	Are directly managed by the program.	

## V Workflow Optimization

### V.1 DICOM Worklist tag mapping

SIDICOM WLS and SIDEXIS communicate with each other via the Plug-in interface.

Here the desired fields from the selected order of the worklist are mapped to the SIDEXIS variables:

```
PermanentDentist=<(0032,1032)> ;Requested Physician  
ImageReason=<(0040,0100)(0040,0007)> ;Scheduled Procedure Step Description  
ImageType=<(0040,0100)(0040,0008)(0008,0100)> ;Scheduled Action Item Code Value
```

These fields can be retrieved from any tag of the worklist order (SPS).

### V.2 Action Item Mapping

#### V.2.1 What does "Action Item Mapping" mean and when is it useful?

Action Item Mapping allows to "translate" customer-specific order codes or common predefined image sequence requests (e.g. bite wing exposures, full mouth series or special panoramic techniques) into a SIDEXIS order request which is customized to the specific capabilities of the Sirona X-ray equipment.

Thus Action Item Mapping means not only specifying orders more accurate – it is the gate to use the full range of the X-ray unit capabilities and to avoid misunderstandings leading to unwanted/unhelpful image types.

In the following the SIDEXIS order request is called "Image Type".

Action Item Mapping takes place in the file sidicom.ini in section [ActionItemMap] and uses the following tags by default:

```
ImageType=<(0040,0100)(0040,0008)(0008,0100)> ;Scheduled Action Item Code Value
```

Dummy values are placed in the ActionItemMap delivered to make input easier. - Excerpt from the "ActionItemMap":

#### [ActionItemMap]

```
IO=XI  
XYZ01=11XI  
XYZ02=12XI  
...  
PAN=XP  
PAN01=01XP  
PAN02=02XP
```

For example, one hospital calls an exposure of the first front tooth in the first quadrant "XYZ0", this name is mapped in SIDEXIS as image type "11XI".

Thus the entries on the right side give all image types that are known in SIDEXIS ("01XP" e.g. stands for the 1st Program of the ORTHOPHOS DS Panoramic device).

The entries on the left side describe the hospital's name for an exposure in the program. This is why it is advisable to replace, wherever possible, the left-hand parts of the "ActionItemMap" section by your own names.

To clarify the best workflow and avoid mistakes we provide in a next chapter some **configuration sheets to be able to discuss the necessary Action Item Mapping with both technical and medical (dental) personnel.**

### V.2.2 Additional comments

In addition to specifying the exact image type **additional comments useful for the X-ray operating personnel** can be placed in an official DICOM tag while producing a DICOM Worklist order.

E.g. tag (0040,0100) Image Reason or (0040,0007) Scheduled Procedure Step can be used.

As our experience shows that in practice mostly tag (0040,0007) Scheduled Procedure Step Description is filled in with such comments we predefined the Image Reason displayed in SIDEXIS as:

ImageReason=<(0040,0100)(0040,0007)> ;Scheduled Procedure Step Description

### V.2.3 Hints for intraoral images

If you want to place orders for more than one single tooth:

Per incoming order in SIDEXIS only one image type/one exposure can be generated.

For creating a FMX-series it is necessary that the Worklist provider already produces image requests for every single image (E.g. for an FMX 14 series one would need 14 single orders in SIDEXIS).

For bite wing exposures we recommend to use the middlest tooth for specifying the order and place comment "Bite wing" into tags (0040,0100) or (0040,0007). This comment is visible in the SIDEXIS order window.

If the tooth number cannot be specified or if this information is not transferred by the worklist please always use "00XI" in Action Item Mapping.

Alternatively you can delete mapping in "ImageType" of the "SLIDA" section and configure SIDEXIS in such a way that the **image type must be determined prior to the IO exposure within SIDEXIS.**

### V.2.4 Action Item Mapping configuration sheets

We provided configuration sheets on the following pages for:

- Intraoral X-rays
- ORTHOPHOS DS
- ORTHOPHOS XG
- Cephalometric X-rays for both ORTHOPHOS DS and XG
- Transversal Slices for both ORTHOPHOS DS and XG

Action Item Mapping configuration sheet for Intraoral images		
Customer-specific X-ray order code	SIDEXIS image type	Shall a DICOM tag be used to specify image reason?
Intraoral-11	11XI	not necessary
Intraoral-12	12XI	not necessary
Intraoral-13	13XI	not necessary
Intraoral-14	14XI	not necessary
Intraoral-15	15XI	not necessary
Intraoral-16	16XI	not necessary
Intraoral-17	17XI	not necessary
Intraoral-18	18XI	not necessary
Intraoral-21	21XI	not necessary
Intraoral-22	22XI	not necessary
Intraoral-23	23XI	not necessary
Intraoral-24	24XI	not necessary
Intraoral-25	25XI	not necessary
Intraoral-26	26XI	not necessary
Intraoral-27	27XI	not necessary
Intraoral-28	28XI	not necessary
Intraoral-31	31XI	not necessary
Intraoral-32	32XI	not necessary
Intraoral-33	33XI	not necessary
Intraoral-34	34XI	not necessary
Intraoral-35	35XI	not necessary
Intraoral-36	36XI	not necessary
Intraoral-37	37XI	not necessary
Intraoral-38	38XI	not necessary
Intraoral-41	41XI	not necessary
Intraoral-42	42XI	not necessary
Intraoral-43	43XI	not necessary
Intraoral-44	44XI	not necessary
Intraoral-45	45XI	not necessary
Intraoral-46	46XI	not necessary
Intraoral-47	47XI	not necessary
Intraoral-48	48XI	not necessary
Bite wing	00XI	Yes, specify the exact region or simply "Bite wing"
Full Mouth series ("FMX")	00XI	Yes, specify the exact series description or common short, e.g. "FMX14"
<b>Notes:</b>		
The above items are image type descriptions to clarify the image types in general. The descriptions can be replaced by the corresponding custom specific image codes. We recommend to enter the codes for every Worklist order created into DICOM tag (0008,0100); Scheduled Action Item Code Value.	This code has to be used in "SIDICOM.ini" to map it to the customer-specific code. For SIDEXIS versions lower than XG this code appears also in the SIDEXIS GUI but w/o any verbal description.	We recommend to use DICOM tag (0040,0007) to further specify the image reason

**Action Item Mapping configuration sheet for ORTHOPHOS DS, panoramic part**

Customer-specific X-ray order code	SIDEXIS image type	Program name/icon SIDEXIS	Shall a DICOM tag be used to specify image reason?
Standard panoramic exposure	01XP	Normal OPG (P1)	Not necessary
Standard panoramic exposure, restricted to teeth, w/o ascending branches	02XP	Normal OPG, reduced (P2)	Not necessary
Sinus maxillaris, 2-on-1 film subdivision	03XP	Paranasal sinuses (S2)	Not necessary
Lateral TMJ with ascending branches	04XP	Temporomandibular joints lateral, ascending branches (TM3)	Not necessary
TMJ in posterior/anterior view	05XP	Temporomandibular joints p.a., ascending branches (TM4)	Not necessary
TMJ with mouth open and closed, on one image	06XP	Lateral view of the temporomandibular joints, closed/open mouth (TM1)	Not necessary
TMJ in posterior/anterior view, with mouth open and closed on one image	07XP	P.a. view of the temporomandibular joints, closed/open mouth (TM2)	Not necessary
Multi-layer exposures of the lateral TMJ	08XP	Lateral multi-layer image of the temporomandibular joints (TM5)	Not necessary
Multi-layer exposures TMJ in posterior/anterior view	09XP	P.a. multi-layer image of the temporomandibular joints (TM6)	Not necessary
Standard pediatric panoramic exposure (with considerable dose reduction)	10XP	Normal OPG for children (P10)	Not necessary
Panoramic exposure with constant 1,25x magnification (e.g. for implantology)	11XP	Normal OPG with constant magnification (P1)	Not necessary
Panoramic exposure with increased layer thickness in the anterior region	12XP	Thick layer, anterior region (P12)	Not necessary
Paranasal tomographic exposures e.g. blow-out fractures	13XP	Maxillary sinuses, 2 images (S1)	Not necessary
Standard panoramic exp., left half side for dose reduction in control exposures	14XP	Normal OPG Left (P1)	Not necessary
Standard panoramic exp., right half side for dose reduction in control exposures	15XP	Normal OPG Right (P1)	Not necessary
Multi-layer in the molar teeth area	16XP	Multi-layer image, lateral teeth (TS1)	Not necessary
<b>Notes:</b>			
<p>The items above are image type descriptions to clarify the image types in general. The descriptions can be replaced by the corresponding custom specific image codes.</p> <p>We recommend to enter the codes for every created Worklist order into DICOM tag (0008,0100);Scheduled Action Item Code Value.</p>	<p>This code has to be used in "SiDICOM.ini" to map it to the customer-specific code. For SIDEXIS versions lower than XG this code appears also in the SIDEXIS GUI w/o any verbal description.</p>	<p>From SIDEXIS XG version 1.3x onwards there are no program codes for ORTHOPHOS DS/Ceph visible. All codes appearing in brackets are valid for the new ORTHOPHOS XG generation only. The customer can use the highlighted program icons in the lower range of the image type dialog to identify the program codes. These icons are the same as labelled on the panel of the unit if display of ORTHOPHOS XG icons is switched off in SIDEXIS.</p>	<p>No DICOM tag needs to be used to further specify the image reason.</p>

**Action Item Mapping configuration sheet for ORTHOPHOS XG, panoramic part**

Customer-specific X-ray order code	SIDEXIS image type	Program name/icon SIDEXIS	Shall a DICOM tag be used to specify image reason?
Standard panoramic exposure	01XP	Normal OPG (P1)	Yes, please specify if you want artefact-free or constant magnification
Standard panoramic exposure, restricted to teeth, w/o ascending branches	02XP	Normal OPG, reduced (P2)	Yes, please specify if you want artefact-free or constant magnification
Sinus maxillaris, 2-on-1 film subdivision	03XP	Paranasal sinuses (S2)	Not necessary
Lateral TMJ with ascending branches	04XP	Temporomandibular joints lateral, ascending branches (TM3)	Not necessary
TMJ in posterior/anterior view	05XP	Temporomandibular joints p.a., ascending branches (TM4)	Not necessary
TMJ with mouth open and closed, on one image	06XP	Lateral view of the temporomandibular joints, closed/open mouth (TM1)	Not necessary
TMJ in posterior/anterior view, with mouth open and closed on one image	07XP	P.a. view of the temporomandibular joints, closed/open mouth (TM2)	Not necessary
Multi-layer exposures of the lateral TMJ	08XP	Lateral multi-layer image of the temporomandibular joints (TM5)	Not necessary
Multi-layer exposures TMJ in posterior/anterior view	09XP	P.a. multi-layer image of the temporomandibular joints (TM6)	Not necessary
Standard pediatric panoramic exposure (with considerable dose reduction)	10XP	Normal OPG for children (P10)	Yes, please specify if you want artefact-reduced or constant magnification
Panoramic exposure with constant 1,25x magnification (e.g. for implantology)	11XP	Normal OPG with constant magnification (P1)	Yes, please specify that you want constant magnification
Panoramic exposure with increased layer thickness in the anterior region	12XP	Thick layer, anterior region (P12)	Not necessary
Paranasal tomographic exposures e.g. blow-out fractures	13XP	Maxillary sinuses, 2 images (S1)	Not necessary
Standard panoramic exp., left half side for dose reduction in control exposures	14XP	Normal OPG Left (P1)	Yes, please specify that you want restrict it to left half side
Standard panoramic exp., right half side for dose reduction in control exposures	15XP	Normal OPG Right (P1)	Yes, please specify that you want restrict it to right half side
Multi-layer in the molar teeth area	16XP	Multi-layer image, lateral teeth (TS1)	Yes, please specify the molar of interest
<b>Notes:</b>			
The above items are image type descriptions to clarify the image types in general. The descriptions can be replaced by the corresponding customer-specific image codes. We recommend to enter the codes for every created Worklist order into DICOM tag (0008,0100); Scheduled Action Item Code Value	This code has to be used in "SiDICOM.ini" to map it to the customer-specific code. For SIDEXIS versions lower than XG this code appears also in the SIDEXIS GUI w/o any verbal description.	Since SIDEXIS XG version 1.3x all codes appearing in brackets are representing the main program to be selected at the Easypad for the new ORTHOPHOS XG Generation. Additional information like artefact-free, constant magnification etc. must be derived from field "Image reason" when orders are accepted inside SIDEXIS.	We recommend to use DICOM tag (0040,0007) to further specify the image reason.

## Action Item Mapping configuration sheet – Cephalometric X-rays for both ORTHOPHOS DS and XG

Customer-specific X-ray order code	SIDEXIS image type	Program name/icon SIDEXIS	Shall DICOM tag (0040,0007) be used to specify image reason?	ORTHOPHOS XG program
Skull exposure p/a	01XC	P.a. view of skull	Not necessary	C1
Skull exposure a/p	02XC	A.p. view of skull	Not necessary	C2
Lateral skull exposure	03XC	Lateral view of skull	Not necessary	C3
Carpus exposure	04XC	Hand	Not necessary	C4
<b>Notes:</b>				
The above items are image type descriptions to clarify the image types in general. The descriptions can be replaced by the corresponding customer-specific image codes. We recommend to enter the codes for every created Worklist order into DICOM tag (0008,0100); Scheduled Action Item Code Value.	This code has to be used in "SiDICOM.ini" to map it to the customer-specific code. For SIDEXIS versions lower than XG this code appears also in the SIDEXIS GUI w/o any verbal description.			

**Action Item Mapping configuration sheet Transversal Slices for both ORTHOPHOS DS and XG**

Customer-specific X-ray order code	SIDEXIS image type	ORTHOPHOS program	SIDEXIS program	Shall a DICOM tag be used to specify image reason?
Transversal Slice right molars upper jaw, teeth 18/17/16	00XS or ??XS	P17	TSA UJ Molars 18/17/16	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left molars upper jaw, teeth 28/27/26	00XS or ??XS	P18	TSA UJ Molars 28/27/26	Yes, if custom specific order code does not contain the teeth number
Transversal Slice right molars lower jaw, teeth 48/47/46	00XS or ??XS	P19	TSA LJ Molars 48/47/46	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left molars lower jaw, teeth 38/37/36	00XS or ??XS	P20	TSA LJ Molars 38/37/36	Yes, if custom specific order code does not contain the teeth number
Transversal Slice right premolars upper jaw, teeth 15/14	00XS or ??XS	P21	TSA UJ Premolars 15/14	Yes, if custom specific order code does not contain the teeth number
Transversal Slice right canine upper jaw, tooth 13	00XS or 13XS	P21	TSA UJ Canine tooth 13	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left premolars upper jaw, teeth 25/24	00XS or ??XS	P22	TSA UJ Premolars 25/24	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left canine upper jaw, tooth 23	00XS or 23XS	P22	TSA UJ Canine tooth 23	Yes, if custom specific order code does not contain the teeth number
Transversal Slice right premolars lower jaw, teeth 45/44	00XS or ??XS	P23	TSA LJ Premolars 45/44	Yes, if custom specific order code does not contain the teeth number
Transversal Slice right canine lower jaw, tooth 43	00XS or 43XS	P23	TSA LJ Canine tooth 43	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left premolars lower jaw, teeth 35/34	00XS or ??XS	P24	TSA LJ Premolars 35/34	Yes, if custom specific order code does not contain the teeth number
Transversal Slice left canine lower jaw, tooth 33	00XS or 33XS	P24	TSA LJ Canine tooth 33	Yes, if custom specific order code does not contain the teeth number
Transversal Slice anterior teeth upper jaw, teeth 12/11/21/22	00XS or ??XS	P25	TSA UJ Front 12/11/21/22	Yes, if custom specific order code does not contain the teeth number
Transversal Slice anterior teeth lower jaw, teeth 32/31/41/42	00XS or ??XS	P26	TSA LJ Front 32/31/41/42	Yes, if custom specific order code does not contain the teeth number
<b>Notes:</b>				
The above items are image type descriptions to clarify the image types in general. The descriptions can be replaced by the corresponding customer-specific image codes. We recommend to enter the codes for every created Worklist order into DICOM tag (0008,0100); Scheduled Action Item Code Value.	Replace "??" by the tooth number. If tooth is not specified replace "??" by "00". e.g. "16XS" or "00XS" if tooth is not specified.	User must select it depending on tooth number	Items valid only for SIDEXIS XG version 1.3x or higher	If you can manage to specify the tooth number for every single tooth it can be taken over directly into the SIDEXIS image type and a DICOM tag needs not to be used. Otherwise we recommend to use DICOM tag (0040,0007); Image Reason.

## VI Operational Hints (Questions & Answers, "Q & A")

**Q** Why are more orders displayed in the worklist view than specified in the query dialog at "Limit query results to ... items"?

**A** It is the result of the answers to the query coming back from the worklist provider (MWL SCP) that is limited rather than the length of the list displayed. Orders somehow being processed may still remain displayed.

**Q** Why is the number of orders displayed in the worklist viewer ever rising up when I repeat a query several times?

**A** Orders already being somehow in processing may still remain displayed as well. Furthermore old query result orders not yet being processed are deleted only when option "DelOldQueryResults" is set in "SiDICOM.ini". How many orders are added by every subsequent query and if it be new ones at all is depending on the usage of the MWL SCP in the meantime by personnel and on the broadness of the query, of course.

**Q** An image was made by mistake and shall not be sent to the DICOM archive. What is to be done?

**A** Assumed that no automatic sending is configured, the Worklist order can be reactivated at any time and a further image can be made which then can be sent manually.

**Q** Some of the images in the DICOM archive contain neither order data nor Accession no.. What went wrong?

**A** Possibly the images are sent with the command "Create New DICOM Study". This command is not based on worklist data. A new Study based on SIDEXIS is generated. The accession number must be filled in manually in a popup dialog. Or an image exposure is made without further transmitting the order to SIDEXIS.

**Q** I found that dialog "Job finished successfully" is not visible sometimes. How do I know if the images were sent properly to DICOM archive?

**A** Sometimes dialog "Job finished successfully" is suppressed depending on the operating system used. This is a problem caused by Microsoft. At least images sent successfully will be marked as "Finished" within the worklist view.

**Q** Sometimes a Popup asks me to enter the Accession No. (or other data). What is to be done?

**A** First enter the data wanted if you know it. Otherwise please ask your technical support for clarification. It may happen that this is due to the operational workflow used or that a purely technical issue occurred.

## VII Technical Support

### VII.1 Logging of SIDICOM WLS

SIDICOM WLS is equipped with a distinct logging interface to search for information and errors.

a)

All logging information is stored in subdirectory ".\SIDICOM\Log". The most important file is "SIDICOM.log". Here all activities and error messages of SIDICOM WLS are recorded. It is thus the first place to consult when facing any problem, e.g. communication or link problems.

b)

If "Verbose" is switched on, additional Logging is made (particularly internal states).

c)

Being in the "Debug" mode, messages of the DICOM communication to be sent and sometimes also messages to be received are additionally stored as text files in the log directory. These might be of help when facing an incorrectly inserted DICOM tag. Depending on the used service the files will be named like, e.g. "C\_Find\_RQ.txt" or "C\_Find\_RSP.txt". In addition there are detailed communication logs in "SIDICOM.log".

d)

A file "Unknown Command.txt" will be created in ".\SIDICOM\Log" in case that the service flag in "SIDICOM.ini" is set to create files or in case debugging is activated. This is just a message but not an error.

e)

#### Settings in "SIDICOM.ini":

[DEBUG]

Debug=0	⇒ Switching "debug" mode on or off, set to 1 only during trouble-shooting.
Verbose=0	⇒ Switching verbose mode on or off, set to 1 only during trouble-shooting.
LogFileOn=1	⇒ SIDICOM logging, default = 1; should not be switched off.
LogFileLen=5000	⇒ Maximum length of the SIDICOM log file, default = 5.000 kB. (for new installations since V.3.0, otherwise 1.000 kB.)

## VII.2 Error messages

The following table gives you a list of error messages of SIDICOM WLS and offers you possible solutions:

Error message	Description	Proposal of solution
Error: Wrong installation. Please check your DLL versions.	The installation of SIDICOM WLS is incorrect.	Save old installation
Error: The DICOM message cannot be sent permanently. Continue trying?	The transfer of images or PPS messages to the DICOM partner repeatedly failed. This happens if the partner is not available at all or for a longer period of time.	This error message comes always together with a detail error message. Eliminate the cause of the detail error. Check the network connection and availability of partner.
Error: Unable to open the association.	The connection with the partner could not be established.	Look for detailed error descriptions in the log files.  Check the physical network connection. Examine the TCP/IP stack. Examine the configuration of the partner. (AE title, IP address, port) in "SIDICOM.ini" file. Try a "ping" and a DICOM "echo test".  Check the access rights of SIDICOM at the DICOM partner.
Error: No presentation context for this SOP class. Please check logfile for detail information.	During the send of image data to the C-Store partner, an error occurs. Possibly, the C-Store partner does not support the image type.	Check logfiles for a detailed error description.  Check the C-Store-partner's image type support in the configuration dialog, respectively in the partner section in "SIDICOM.ini", and change it if necessary.
Error while receiving the worklist.	While receiving worklist data errors occurred.	Check the section [WQuery] in SIDICOM.ini and the query responses of the worklist SCP on DICOM conformance.
Error: Unexpected response.	SIDICOM WLS has received an unknown message from the DICOM partner.	Check the response of the partner for DICOM conformance.
Error while creating a file. Error while creating an empty file. Error while creating an empty file.	SIDICOM WLS cannot create a file. This error can occur when access rights are missing or when using read-only files.	Please check access rights on the paths of the file system which are used by SIDICOM WLS. Remove any read-only protection of the file or the disk. Restart the system.
Other errors	Other errors occurring repeatedly can not be described here in general..	Please contact the Sirona hotline. If an error has occurred for the first time, check your installation and configuration first.

### VII.3 Technical hints (Questions & Answers, "Q & A")

**Q** The **number of orders displayed** in the worklist viewer rising when I repeat a query several times. What is wrong?

**A** Please assure that queries were performed by entering a date where orders are valid. Furthermore old query result orders not yet being processed are deleted only when option "DelOldQueryResults" is set in "SIDICOM.ini". How many orders are added by every subsequent query and if it be new ones at all is depending on the usage of the MWL SCP in the meantime by personnel and on the broadness of the query, of course.

**Q** **Worklist orders** will not be processed properly or images sent do not contain all expected data. What is wrong?

**A** Please find below some hints which may help you to identify the problem:

- We recommend activating automatic sending of images to DICOM for the default SCP partner. If this can not be implemented, please assure that images will be sent to DICOM manually.
- Closing SIDEXIS will cause loss of all unprocessed orders - a new transmit of the order within SIDICOM to SIDEXIS is necessary.
- Images sent manually will contain the order data (like Accession No. etc.) only in case that orders were accepted and processed without intermittently closing of images/application or intermittently switching patients.
- Do not double-click the orders in SiDICOM several times. This produces unwanted X-ray orders in SIDEXIS confusing the user on account of differently expected patient sequence.
- Do not use the "SIDEXIS" button in SIDICOM for task switching – it is an order transmission button.

**Q** How can I read out from the DICOM image **which SIDEXIS image type was performed?**

**A** DICOM tag (0008, 2111) "Derivation Description" contains a text description of how this image was derived, e.g. "01XP". DICOM tag (0018, 1508) contains enumerated values for the positioner Type, e.g. "NONE" for intraoral X-rays or "CEPHALOSTAT" / "PANORAMIC" for ORTHOPHOS X-rays.

**Q** We do not have established a Worklist connection yet but we already want to send e.g. the **Accession No.** to the DICOM archive. Is that possible?

**A** Yes, a window will pop up prompting for a missing Accession No. automatically in standard installation. If you want to enter different data please refer to chapter "Monitoring tag values" and edit section "[CheckImageHeader]" of "SiDICOM.ini" according to manual.

**Q** We have established a Worklist connection yet and we do not need to send the Accession No. to the DICOM archive. But every time **an irritating window pops up**. Is that possible?

**A** Yes, we recommend to disable tag monitoring to avoid the dialog which asks for entering the Accession No.. Please delete or comment out line "(0008,0050);Accession Number" in section "[CheckImageHeader]" of "SiDICOM.ini".

**Q** Please tell me more about **automatic image transfer** to DICOM partners!

**A** Automatic image transfer to the default DICOM partner can be activated in the configuration dialog of SIDICOM WLS.

**Q** How can I **verify the data output** of the Sirona application?

**A** In case you just want to verify single tags we recommend usage of the built-in tag monitoring. Configuration takes place in section "[Check ImageHeader]" of "SIDICOM.ini". A dialog prompting for otherwise empty tags will pop up.

**Q** The actual IP of my computer is unknown to me. May I use the **loopback IP** for verifying connectivity to other applications running on the same computer?

**A** No, the loopback IP does not always work properly in these cases. Please use the NIC 's (Network Interface Card) IP.

**Q** I reviewed the configuration files and wonder what happens after x **GarbageDays**?

**A** The content of directory ".\SIDICOM\SPS\" will be cleaned up from older entries and 'content table' "SPSList.dat" will be modified accordingly. So it is not simply a deletion of older files.

**Q** A technical issue is to be solved. How can I **delete existing orders** in the Sirona application?

**A** Please simply delete the contents of directory ".\SiDICOM\SPS\". A new 'content table' "SPSList.dat" will be generated automatically on next startup. Attention: All information still existing about orders already sent/processed of the worklist view will be lost as well, of course!

**Q** I doubt that the image resolution was transferred correctly into the DICOM archive. How to check?

**A** According to the DICOM Standard resolution information is forwarded differently depending on the image format used.

For the Secondary Capture (SC) format the resolution can be found in tag (0028,0030) Pixel Spacing. For the DX format it is to be found in tag (0018,1164) Imager Pixel Spacing. Please verify that your viewer etc. is recognizing tag (0018,1164) properly.

**Q** I set a worklist order to "Finished" or "Discontinued" manually. Why does the PPS partner not receive a message?

**A** Manually setting worklist orders is basically an orientation for the user in the worklist. Notifying the PPS SCP in this case is only possible if an X-ray exposure has already been made and been transmitted from SIDEXIS. This may be useful if there was a sending error in the past and the image has already been transmitted to the archive through other ways.

## VII.4 Support Contact

Support is provided by our hotline:

Tel. +49 (0) 62 51 / 16 – 16 16

Fax: +49 (0) 62 51 / 16 – 18 18

Email: [product.service@sirona.de](mailto:product.service@sirona.de)

If possible, keep at hand the logging and configuration files as well as information on your network environment and your communication partner for an easier troubleshooting.

**Providing the complete content of the local "C:\Program Files\Sirona\SIDICOM WLS" respectively "C:\SiDICOM" directory and all "\*.ini" files from the SIDEXIS directory would be of great advantage.**

At least the following files shall be provided:  
All files with extension ".txt", ".log", ".ini" from the installation directory and subdirectories and from the SIDEXIS directory!

---

We reserve the right to make any alterations which may be due to technical improvements.

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© Sirona Dental Systems GmbH 2002 - 2007
Sirona Dental Systems LLC
**Sirona Dental Systems GmbH**